austria

portogallo

rumeno

turco

arabo

GB

Medical Device Class II b  
INTERNATIONAL PATENT

## GENGIGEL® forte Oral Gel

*Ideal for application to wounds and during the post-operative phase.  
Specifically designed to ensure simple, rapid, efficient use.*

Hyaluronic acid is the most physiologically important proteoglycan found naturally in the mucosal extracellular matrix, where it interacts with other components of connective tissue to provide the strength and elasticity essential for good oral health.

Hyaluronic acid has a number of beneficial properties that contribute to a successful programme of dental treatment and oral hygiene.

In particular, hyaluronic acid can absorb 50 times its own dry weight of water.

This characteristic enables hyaluronic acid to act as a barrier to the penetration of micro-organisms and other toxins which are likely to contribute to the aetiology of periodontal disease.

Gengigel® Forte

- Contains a non-animal derived, biotechnologically manufactured source of high quality, high purity hyaluronic acid in the form of a high molecular weight sodium salt.
- Has a long-lasting local action.
- Adheres to the mucous membrane due to the special nature of the product formulation.
- Can be used, after medical advice, without restriction by all patients including children, pregnant women, the elderly and diabetics.

### Indications

Gengigel® Forte

- To be used on wounds that can heal by secondary intent.
- Reduces pain and discomfort, especially after traumas, surgical or periodontal work, implants, tooth extractions and subgingival root planing.
- Promotes and accelerates the anti-oedematous and anti-inflammatory processes.

### Directions For Use

Immediately after dental treatment apply a small quantity of gel over the surface of the affected area or, after washing hands, thoroughly apply with a finger and massage it gently on to the gum until the affected area is entirely covered.

Use 1-2 times a day, or as required, unless otherwise advised by a doctor.

### Composition

Principal ingredient: Sodium Hyaluronate

Other components: Aqua, Xylitol, Cellulose Gum, Alcohol, PEG 40 Hydrogenated Castor Oil, Glycerophosphocholine, Sodium Phosphate, Polyvinyl Alcohol, Carbomer, Dichlorobenzyl Alcohol, Trisodium phosphate, Sodium Lactate, Sodium Hydroxide, Citric Acid.

### Packaging

Gengigel® Forte oral gel tube 8 ml

**Contraindications:** the use of GENGIGEL® Forte is contraindicated in the event of hypersensitivity to any of its ingredients.

**Precautions and warnings:** do not use in case of any visible damage. Store at room temperature, protect against direct heat sources. Do not use after the expiry date. Keep out of children's reach. Should the symptoms persist, consult a healthcare professional.

A

Medical Device Class II b  
INTERNATIONAL PATENT

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**(SLO)** **Medical Device Class II b**  
INTERNATIONAL PATENT

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**Sestavine**

Glavna sestavina: Sodium Hyaluronate

Ostale sestavine: Aqua, Xylitol, Cellulose Gum, Alcohol, PEG 40 Hydrogenated Castor Oil, Glycerophosphocholine, Sodium Phosphate, Polyvinyl Alcohol, Carbomer, Dichlorobenzyl Alcohol, Trisodium phosphate, Sodium Lactate, Sodium Hydroxide, Citric Acid.

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**(SRB)** **Medical Device Class II b**  
INTERNATIONAL PATENT

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**Sastojci**

Glavni sastojak: Sodium Hyaluronate

Ostali sastojci: Aqua, Xylitol, Cellulose Gum, Alcohol, PEG 40 Hydrogenated Castor Oil, Glycerophosphocholine, Sodium Phosphate, Polyvinyl Alcohol, Carbomer, Dichlorobenzyl Alcohol, Trisodium phosphate, Sodium Lactate, Sodium Hydroxide, Citric Acid.

**Pakovanje**

Gengigel® **Forte** peroralni gel: tuba od 8 ml

**Kontraindikacije**

Primena se ne preporučuje kod poznate preosetljivosti na bilo koji sastojak.

**Mere opreza i upozorenja**

Ne upotrebljavajte proizvod ako na njemu primetite vidljiva oštećenja. Čuvajte na sobnoj temperaturi, zaštićeno od izvora toplote. Ne upotrebljavajte posle datuma isteka roka upotrebe. Čuvajte van domašaja dece. Ukoliko se simptomi ne povuku, posavetujte se sa stomatologom.

Broj i datum rešenja o upisu u registar u RS:

**(HR)** **Medical Device Class II b**  
INTERNATIONAL PATENT

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**(H)** **Medical Device Class II b**  
INTERNATIONAL PATENT

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